



IV Flow Controller - Controflow

Intended Use

I. V. Flow Controller is used for adjusting infusion flow speed in-vitro. I. V. Flow Controller is needed for controlled flow of infusion fluids of pediatric patients, critical drug infusions, and during anaesthesia.

Components of The Product

Medbar brand Controflow consists of:

- Hose (tube)
- Lid
- Luer Locks
- Y-Port
- Scale

Product Features

- Gravity type flow controller.
- Constant flow rate.
- Universal male and female luer connections.
- Needle & Needle-free Y-Port injection options available.
- Different tubing length available.

Shelf Life

5 years.

Storage Conditions

Must be stored in dry, ventilated environments, and away from sunlight.



LATEX - FREE

STERILE EO